

# Diabetic Prediction Framework: A Privacy-Enhanced Approach with Metaheuristic Optimization and Blockchain

## 1. Project Overview

This project develops a cutting-edge **Diabetic Prediction System** that integrates advanced machine learning techniques with blockchain technology and a user-friendly web interface. Its primary goal is to provide a robust, accurate, and transparent solution for assessing diabetes risk based on various patient health parameters.

The system leverages an **ensemble** of powerful machine learning models, enhanced by **metaheuristic optimization** for feature selection, to classify individuals into Non-Diabetic, Pre-Diabetic, or Diabetic categories. Crucially, every prediction made is securely recorded onto a custom-built **blockchain**, ensuring an immutable and auditable history of all health assessments. The entire solution is made accessible through an intuitive Flask-based web interface, facilitated by Ngrok for public access, allowing users to input health parameters and receive real-time predictions with verifiable integrity.

## 2. Key Features

- **Accurate Diabetes Prediction:** Utilizes an ensemble of machine learning models (XGBoost and Logistic Regression) for high predictive accuracy.
- **Optimized Feature Selection:** Incorporates **Metaheuristic Optimization** (specifically PSO - Particle Swarm Optimization) to select the most relevant features, enhancing model performance and efficiency.
- **Class Imbalance Handling:** Employs SMOTE (Synthetic Minority Over-sampling Technique) during training to address skewed class distributions, ensuring fair predictions across all diabetes stages.
- **Immutable Prediction Records:** Integrates a custom **Blockchain** to securely store every prediction, providing a tamper-proof and auditable history.
- **User-Friendly Web Interface:** A Flask-based web application with a clean HTML/Tailwind CSS frontend for easy patient data input and result visualization.
- **Real-time Risk Assessment:** Provides instant predictions along with confidence scores and tailored health recommendations.
- **Public Accessibility:** Deployed using **Ngrok** for easy access and demonstration from any web browser.

## 3. Technologies Used

- **Frontend:**
  - HTML
  - Tailwind CSS
  - Jinja Templating
- **Backend:**
  - Python (Core language for all logic)
  - Scikit-learn (Machine Learning models, preprocessing)
  - NumPy (Numerical operations)
  - Pandas (Data manipulation)
  - Joblib (Model persistence)
  - Custom Python Blockchain Implementation (Block and Blockchain classes)
- **Frameworks:**
  - Flask (Web application framework)
  - Pyngrok (For creating public tunnels)

## 4. Project Structure

Diabetic\_Prediction\_Model/

```

├── app.py           # Main Flask application logic
├── blockchain.py    # Custom Blockchain implementation
├── diabetic.csv.csv # Original dataset (ensure this is in your Drive)
├── ensemble_model.pkl # Trained ensemble ML model
├── gender_encoder.pkl # LabelEncoder for 'Gender'
├── original_feature_names.pkl # List of feature names used by the model
├── scaler.pkl       # StandardScaler for feature scaling
├── selected_columns.pkl # Indices of features selected by PSO
├── diabetes_blockchain.json # Persistent storage for the blockchain
├── templates/      # HTML templates for the web interface
│   ├── index.html  # Main prediction input form and results display
│   └── history.html # Blockchain history viewer

```

*(Note: Your Colab notebook file, e.g., YourProjectNotebook.ipynb, would typically reside outside this Diabetic\_Prediction\_Model folder, often in Colab Notebooks/.)*

## 5. How to Run the Project (in Google Colab)

1. **Upload diabetic.csv.csv:** Ensure your diabetic.csv.csv dataset is uploaded to your Google Drive, preferably in the same Diabetic\_Prediction\_Model folder or accessible via '/content/drive/MyDrive/Colab Notebooks/Diabetic/diabetic.csv'.
2. **Open Your Colab Notebook:** Open your project's .ipynb file in Google Colab.

3. **Run Cells Sequentially:** Execute all code cells in your Colab notebook from top to bottom.
  - **Cell 1:** Creates templates directory and writes blockchain.py.
  - **Cell 2:** Writes templates/index.html.
  - **Cell 3:** Writes templates/history.html.
  - **Cell 4:** Writes app.py (contains the ML logic, blockchain integration, and Flask routes).
  - **Cell 5:** Installs libraries, mounts Google Drive, copies all necessary files (app.py, blockchain.py, templates/) to your Diabetic\_Prediction\_Model folder on Drive, changes the working directory, sets up Ngrok, and launches the Flask application.
    - **Important:** You will need to replace "YOUR\_NGROK\_AUTH\_TOKEN\_HERE" with your actual Ngrok authentication token in Cell 5.
    - The Flask app will run on http://localhost:5000 and Ngrok will expose it publicly.
4. **Access the Web App:** Once Cell 5 finishes executing, it will print a Flask App Public URL: https://... link. Click this link to open your diabetes prediction web application in your browser.

## 6. Sample Usage

Navigate to the public URL provided by Ngrok. You can input patient data into the form fields.

### Example Inputs:

- **Non-Diabetic (N):**
  - Gender: F, AGE: 30, Urea: 3.0, Cr: 55, HbA1c: 5.0, Chol: 4.0, TG: 1.0, HDL: 1.5, LDL: 2.0, VLDL: 0.5, BMI: 22.0
- **Pre-Diabetic (P):**
  - Gender: M, AGE: 45, Urea: 4.5, Cr: 70, HbA1c: 6.0, Chol: 5.0, TG: 1.8, HDL: 1.0, LDL: 2.8, VLDL: 0.8, BMI: 26.0
- **Diabetic (Y):**
  - Gender: F, AGE: 60, Urea: 7.0, Cr: 90, HbA1c: 7.5, Chol: 6.0, TG: 3.0, HDL: 0.8, LDL: 3.5, VLDL: 1.2, BMI: 32.0

After submitting, the prediction, confidence, and recommendations will be displayed. You can then click "View History" to see the immutable blockchain record of your prediction.

## 7. Resetting Blockchain History

To clear the `diabetes_blockchain.json` file and reset the prediction history, you can run the following command in a new Colab cell (after mounting Drive and setting `PROJECT_DIR_ON_DRIVE`):

```
import os
BLOCKCHAIN_FILE = os.path.join(PROJECT_DIR_ON_DRIVE,
'diabetes_blockchain.json')
if os.path.exists(BLOCKCHAIN_FILE):
    os.remove(BLOCKCHAIN_FILE)
    print(f"Removed existing blockchain file: {BLOCKCHAIN_FILE}")
else:
    print(f"Blockchain file not found at {BLOCKCHAIN_FILE}. Nothing to remove.")
```

After running this, restart your Colab runtime and re-run all cells to start with a fresh blockchain.

## 8. Future Enhancements

- Implement user authentication and role-based access control.
- Explore more advanced deep learning models for prediction.
- Integrate with a decentralized storage solution (e.g., IPFS) for patient data, linking hashes to the blockchain.
- Enhance the web interface with interactive visualizations of patient data and prediction trends.
- Implement a peer-to-peer network for the blockchain to achieve true decentralization.